TrHLINE

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EN

TrHLINE

HIGH MOLECULAR WEIGHT HYALURONIC ACID SOLUTION

DESCRIPTION

TrHLINE is a viscoelastic solution of high molecular weight (1.9-2.2 MDa) sodium hyalurate obtained from biofermentation source, highly purified and no inflammatory, transparent, sterile, free of endotoxins dissolved in an isotonic and physiologic buffer.

PACKAGING

TrHLINE is supplied in a sterile, non-pyrogenic disposable prefilled Luer-lock syringe containing 2 mL of solution.

INTENDED USE

TrHLINE is intended to be injected intra-articularly in the synovial space of the knee as viscosupplementation medium for synovial fluid when its viscosity has been reduced due to osteoarthritis process.

CONTRAINDICATIONS

Do not administer to patients with known hypersensitivity (allergy) to hyaluronate preparations and in the case of infections or skin diseases in the area of the injection site. The use of **TrHLINE** in any other joints besides the knee and for other conditions such as radiographic OA or any non-OA causes of joint pain is not indicated.

PROPERTIES

TrHLINE provides a lubricant effect on osteoarthritic synovial joints in vivo. **TrHLINE** decreases the pain and increases the functionality of the affected joint. **TrHLINE** is reabsorbed by organism in vivo.

COMPOSITION

TrHLINE is formulated with highly purified sodium hyaluronate, pharmacopeia quality salts and water for injection (WFI). Its composition varies in function of the model:

Model:	TrHLINE 1%	TrHLINE 2%
Content	(mg/mL)	(mg/mL)
Sodium hyaluronate	10	20
Disodium phosphate dodecahydrate	0.6	0.6
Sodium dihydrogen phosphate dihydrate	0.05	0.05
Sodium chloride	8	8
Water For Injection	q.s	q.s

INTERACTIONS

Do not use with other products containing quaternary ammonium salts because sodium hyaluronate can precipitate in their presence. In order to prevent any possible interaction avoid the simultaneous administration with other intra-articular products.

ADVERSE EVENTS

Joint effusion, arthralgia, redness, joint warmth, swelling and/or pain may occur sporadically at the site of injection. These are the most common side effects, which are transient and clinically manageable with short term use of ice and NSAIDs. As for any intra-articular treatment, septic arthritis or infection may rarely occur when general precautions for injections are not observed or the site of injection is not aseptic.

WARNINGS AND PRECAUTIONS

TrHLINE has not been tested in pregnant women: special care should be paid if the product is used in pregnant or lactating women. It is recommended to discontinue breastfeeding during treatment. Use in the elderly: the product must be administered with caution, since physiological functions are usually reduced in these patients. Pediatric use: since the safety of this product has not been established in children, special care should be paid if the product is used in children. Do not inject in patients with bleeding disorders, liver diseases and/or venous or lymphatic stasis. The product should not be used if the package is damaged or broken. The product should not be used after the expiration date. The product must be injected strictly into the intra-articular space otherwise it may cause pain. Do not inject into a blood vessel or surrounding tissues. Do not inject in inflamed joints, the product must be administered after the local symptom has been alleviated. To avoid pain after the injection, it should be recommended to keep the injection site in relative rest for 48 hours after the injection. It is recommended that the patient avoids any strenuous physical activities or prolonged (i. e., more than one hour) weight-bearing activities within 48 hours following the intra-articular injection. **STRICT ASEPTIC PROCEDURE SHOULD BE FOLLOWED. SINGLE USE PRODUCT. DO NOT RE-STERILISE**, since re-sterilization may affect the physico-chemical properties of the product and compromise its efficacy and security. **KEEP OUT OF THE REACH OF CHILDREN**. CE marked and Luer-lock fittings accessories (needles) must be used.


DOSAGE AND ADMINISTRATION

Bring to room temperature before injection. **TrHLINE** has to be injected using a suitable sterile needle (21 to 23 G needle is recommended) in the affected joint at weekly intervals for 3 weeks. Further injections may be administered if necessary. The efficacy of the treatment has been observed to last for at least 6 months. The efficacy of a single treatment cycle of less than the indicated 3 injection has not been established, and therefore is not indicated. This product should only be used by qualified personnel, following all the rules regarding aseptic procedure and the injection technique. Remove any joint effusion, if present, before injection. More than one joint may be treated at the same time. Discard any unused portion of the syringe.

DIRECTIONS FOR USE

- Peel off the sterile barrier and remove the syringe from inside.
- Remove the tip cap of the syringe without touching the syringe tip, in order to maintain it sterile.
- Attach a suitable (recommended 21 to 23G) sterile needle to the tip of the syringe. Verify that the needle is properly locked to the Luer-lock adapter of the syringe before using it.
- Remove air from the syringe, if necessary.
- Inject intra-articularly into the synovial space of affected knee, following aseptic injection procedures.
- A total dose of 2 mL of the pre-filled syringe have to be injected in the affected knee once a week, during an interval of 3 weeks, with a total of 3 injections.

RHEOLOGICAL PROFILE

STORAGE

Protect from light and avoid freezing. Store between 5°C and 25°C.

SELL-BY-DATE

Do not use after the date marked on the package.

MANUFACTURED BY

i+Med S.Coop
Parque Tecnológico de Álava
Calle de Albert Einstein, 15. P. 15
01510 Vitoria-Gasteiz, (Álava - Spain)
www.imasmed.com

Segregation and waste management according to current legislation

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